Mask

Full Face Mask (F5AS)

User Manual

Full Face Mask (F5AS)

Thank you for choosing BMC's Full Face Mask (F5AS). The mask is comfortable, flexible and quiet. It is designed to minimize contact with your face, thus ensuring that you feel comfortable during therapy. This user manual provides you with the information you need for the correct use of your mask. The Mask is Not made with natural rubber latex.

Intended Use

The Full Face Mask (F5AS) channels airflow non-invasively to a patient from a positive airway pressure device such as a continuous positive airway pressure (CPAP) or bi-level system.

The Full Face Mask (F5AS) is:

- to be used by adult patients (> 66 lb / 30 kg) for whom positive airway pressure therapy has been prescribed.
- intended for single-patient re-use in the home environment and multi-patient re-use in the hospital / institutional environment.

CAUTION: In the US, Federal law restricts the sale of this device to sale by or on the order of a physician. Please read the entire user manual before use.

Contraindications

This mask should be used with caution for patients with following conditions:


- Facial deformity.
- Injured facial skin which will contact with the mask during use.

Before Using the Mask

⚠ WARNINGS

- The vent holes must be kept clear.
- Only use F5AS Filter spare parts.
- Not applicable for silicone-allergy user.
- The F5AS full face CPAP masks must remove Filter assemblywhen are intended to be used simultaneously with nebulizer medications that are in the air path of the mask / tube.
- Ensure the filter is completely dry before each use.
- The Filter assembly is for single patient use. Replace between patients.
- The F5AS mask can be used without the Filterassembly.
- This mask should only be used with CPAP or bilevel devices recommended by a physician or respiratory therapist.
- **Explanation of Warning:** CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air will be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation. This warning applies to most models of CPAP systems.
- At low CPAP pressures the flow through the exhaled port may be inadequate to clear all exhaled gas from the tube. Some rebreathing may occur.
- To minimize the risk of vomiting during sleep, the patient should avoid eating or drinking for three hours before using the mask. This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
- Do not use the mask without the Non-Rebreathing Valve in place.
- This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.
- Follow all precautions when using supplemental oxygen.
- Oxygen flow must be turned off when the flow generator is not operating, so that unused oxygen does not accumulate within the flow generator enclosure and create a risk of fire.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration varies, depending on the pressure settings, patient breathing pattern, mask, point of application and leak rate.
- The technical specifications of the mask are provided for your clinician to check if it is compatible with the flow generator. If it is used beyond technical specifications or used with incompatible devices, the seal and comfort of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect the function of the flow generator.
- Stop using the Full Face Mask (F5AS) and consult your physician or sleep therapist, if you have ANY adverse reaction to the use of the mask.
- Remove all packaging before using the mask.
- Images shown here are indicative only. If there is inconsistency between the image and actual product, the actual product shall govern.


Components of the Full Face Mask (F5AS)

Note: The Frame assembly is composed of Frame, Panel, Soft-sleeve, Connector and Elbow assembly. The Elbow assembly is composed of Elbow and Valve. The Filter assembly is composed of Filter and Filter cover.

Components of Full Face Mask (F5AS)


Fitting the Mask

Disassembling the Mask

Notes:

1. The elbow assembly, the panel and the connector cannot be dissembled from the frame assembly.
2. The elbow assembly cannot be disintegrated.
3. The valve cannot be dissembled, or damage may be incurred and the re-installation will become more difficult.

Cleaning the Mask at Home

Notes:

- The mask and headgear can only be cleaned through hand washing.
- The elbow assembly, the panel and the connector cannot be dissembled from the frame assembly for washing.
- The elbow and the valve cannot be dissembled from the elbow assembly for washing.
- The Filter assembly should be washed at least once a week. Dry completely before reuse. Simply replace it with a new one after being used for 30 days or after 8 times of cleaning.
- Gently wash the Filter assembly to avoid damaging it. Do not soak for more than 6 minutes.

⚠ CAUTIONS

- Do not use solutions containing bleach, chlorine, alcohol, aromatics, moisturizers, antibacterial agents, or scented oils to clean any part of the system or air tubing. These solutions may cause damage and reduce the life of the product.
- Exposing any part of the system or tubing to direct sunlight or heat may cause deterioration.
- If any visible deterioration of a component is apparent (cracking, crazing, tears etc), the component should be discarded and replaced.

Daily / After Each Use

- Disassemble the mask components according to the disassembly instructions.
- Thoroughly clean the separated mask components (excluding headgear), by gently rubbing in warm (approx. 30°C) water using mild liquid detergent for 10 min.
- Use a soft bristle brush to clean the vent.
- Rinse all components well with drinking quality water and allow to air dry out of direct sunlight.
- When all components are dry, reassemble according to the reassembly instructions.


Weekly

Handwash the headgear and all components in warm (approx. 30°C) and mild liquid detergent for up to 10 minutes. Rinse the components well and allow them to air dry out of direct sunlight.

Reprocessing the Mask between Patients

This mask should be reprocessed when used between patients. Cleaning and disinfection instructions are available from the BMC website, www.bmc-medical.com. If you do not have internet access, please contact your BMC representative.

Reassembling the Mask

Troubleshooting


Problem	Possible Reason	Possible Solution
Mask won't seal properly or is uncomfortable.	Mask may have been fitted incorrectly.	Carefully follow instructions in "Fitting the Mask" section. Make sure the headgear is not over-tightened.
	Mask size is wrong.	Consult your clinician.
Mask leaks around the face.	The cushion is misplaced on the cushion frame.	Check insertion of the cushion and reinsert correctly according to the instructions in "Reassembling the Mask" section.
Mask is too noisy.	Mask size is wrong.	Consult your clinician.
	Vents are blocked or partially blocked.	Clean the vents according to the instructions in "Cleaning the Mask at Home" section.

Technical Specifications

F5AS EXHAUST FLOW

The Full Face Mask (F5AS) has a unique Filter assembly that is designed to diffuse air.

The mask can be used with or without the Filter assembly.

Pressure (cmH ₂ O)	4	10	12	17	24	30
Flow Rate (L/min)	22	37	41	49	59	67
With Filter	24	39	43	51	61	69
With out Filter						

Dead Space Information	Dead space is the empty volume of the mask up to the swivel. The dead space of the mask varies according to cushion sizes but is less than 225 ml.
Therapy Pressure	4 to 30 cmH ₂ O
Resistance	Drop in Pressure measured (average for 3 sizes) at 50 L/min: 0.15 cmH ₂ O, at 100 L/min: 0.5 cmH ₂ O
Inspiratory and Expiratory Resistance	The inspiratory resistance of the mask (in combination with the Non-Rebreathing Valve) is 1.8 cmH ₂ O at 50 L/min. The expiratory resistance of the mask (in combination with the Non-Rebreathing Valve) is 2.0 cmH ₂ O at 50 L/min.
Anti Asphyxia Valve open-to-atmosph here pressure	≤ 2.8 cmH ₂ O
Anti Asphyxia Valve closed-to-atm osphere pressure	≤ 2.8 cmH ₂ O
Sound	DECLARED DUAL-NUMBER NOISE EMISSION VALUES in accordance with ISO 4871. The A-weighted sound power level of the mask with Filter is 25 dBA, with uncertainty 3 dBA. The A-weighted sound pressure level of the mask with Filter at a distance of 1 m is 17 dBA, with uncertainty 3 dBA. The A-weighted sound power level of the mask without Filter is 28 dBA, with uncertainty 3 dBA. The A-weighted sound pressure level of the mask without Filter at a distance of 1 m is 20 dBA, with uncertainty 3 dBA.
Environmental Conditions	Operating temperature: +5°C to +40°C (41°F to 104°F) Operating humidity: 0% ~ 93% relative humidity non-condensing Storage and transport: -25°C to +70°C (-13°F to 158°F) Storage and transport humidity: 0% ~ 93% relative humidity non-condensing

Storage

Ensure that the mask is thoroughly clean and dry before storing it for any length of time. Store the mask in a dry place out of direct sunlight.

Disposal

The mask does not contain any hazardous substances and may be disposed of with your normal household refuse.

Symbols

System and Packaging

⚠	Caution, consult accompanying documents.
⚠	Indicates a Warning or Caution and alerts you to a possible injury or explains special measures for the safe and effective use of the device.
LOT	Lot number
REF	Catalog number
ⓧ	Not made with natural rubber latex
70°C 158°F	Temperature limitation
90%	Humidity limitation
⚡	Manufacturer
EC REP	Authorized Representative in the European Community

Limited Warranty

It is warranted that the Full Face Mask (F5AS), including Frame Assembly, Cushion, Headgear, Swivel and other accessories shall be free from defects in material and workmanship for a period of ninety (90) days from the date of purchase by the initial consumer.

To exercise the rights under this warranty, contact the local authorized dealers.

The expected service life of Full Face Mask (F5AS) is one year.

MANUFACTURER:
BMC MEDICAL CO., LTD.
Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA

URL: en.bmc-medical.com
E-mail: intl@bmc-medical.com
Tel: +86-10-51663880
Fax: +86-10-51663880 Ext. 810

EU AUTHORISED REPRESENTATIVE:
Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

Manufacturing Site:
BMC (Tianjin) Medical Co., Ltd.
2/F North Area and 3/F, Building No.4, No.1 Xinxing Road, Wujing District, (301700) Tianjin, P.R.China
Tel: +86-22-82939881

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